

Wireless Display

Model FR95X-AUX

Operate your scale from anywhere with this our repeater display. This display is the same size and has the same functions as the SRV95X series scale displays.

The display features a desk mounted large 4.3" color touch screen display that is easy to read. It connects to the main display using peer to peer wifi so no network is required.

Display can sit on a desk or mount flush to a wall.

Comes with SR standard 2 year warranty.



DISPLAY

Full function repeater display

COMPATIBILITY

SRV95X Floor Scales

SRV957IFS In-Floor Scales

WIRELESS

Note: facility does not need wi-fi for this feature to work. Scale comes with peer to peer wi-fi to communicate with the primary display.

POWER SUPPLY

5Vdc 2A medical grade power supply

QR WEIGHT DISPLAY

Simply touching the displayed weight will switch the weight to a QR code for easy scanning into your medical record system. Touch again to switch back.

