

## Mobile Platform Scale Model SRV94X Series

Whether you are a mobile or stationary veterinary practice, SR Mobile Platform Scales offer a durable and reliable option. The SRV945, SRV947, and SRV949 are lightweight scales designed for veterinarians and zoo keepers to obtain weights by taking the scale to the animal.

Constructed of high strength aluminum, the SRV94X Series are as durable as they are functional. The scales' surfaces have an extremely durable polyurea coating that offer grip along with ease of cleaning.

Although lightweight, the SRV94X series offers an impressive 400lb weight capacity. Optionally, the scale can be remotely operated using the SR Scales app on your phone or tablet.

### PLATFORM SIZE

• SRV945	24" x 36"
• SRV947	24" x 48"
• SRV949	24" x 60"

### NON-SLIP PLATFORM

Durable and easy to clean polyurea coating over aluminum with aggregate mixed in to create a surface animals are comfortable walking on.

### BATTERY OPERATED

6 "AA" batteries. Auto Power Down (APD) programmable from 60 - 300 seconds.

### OPTIONAL EQUIPMENT

- SRV945CC carry case for SRV945 only.
- Wireless module required for use with SR Scales App (must be factory installed).



SRV947 Shown



SRV945CC optional carry case for SRV945 only.

# Mobile Platform Scale

## Model SRV94X Series


**SR Scales**<sup>®</sup>  
by SR Instruments, Inc.

Quality. Reliability. Accuracy.


### SYSTEM SPECIFICATIONS

<b>MAX WEIGHT</b>	400 lb / 182 kg
<b>DISPLAY RESOLUTION</b>	.2 lb / .1 kg
<b>BUTTON FUNCTIONS</b>	Zero/Weigh, lb/kg, Hold/Recall
<b>DISPLAY UNITS</b>	Pounds / Kilograms
<b>MODEL SIZE</b>	<b>SRV945:</b> 24" x 36" (dim A) <b>SRV947:</b> 24" x 48" (dim A) <b>SRV949:</b> 24" x 60" (dim A)
<b>HOLD</b>	Freezes the weight on the display
<b>POWER</b>	Six (6) AA Batteries
<b>AUTO ZERO</b>	One button operation
<b>SCALE WEIGHT</b>	<b>SRV945:</b> 27 Pounds <b>SRV947:</b> 36 Pounds <b>SRV949:</b> 45 Pounds
<b>WARRANTY</b>	Two (2) Years

### DIMENSIONS in [cm]



### PACKAGING



### ACCESSORIES AND OPTIONS

DESCRIPTION	PART NUMBER
<b>WIRELESS MODULE</b>	Add -BT Suffix to P/N
<b>SRV945 CARRY CASE</b>	SRV945CC

Please note the wireless models still have the built in display as well. The addition of the wireless will also consume more battery and is not recommended for facilities that take a high volume of weights such as a veterinary clinic.

Specifications subject to change without notice.